

Atlas HotTaq 2x PCR Mix - Data Sheet

Cat. No.	Pack Size
BAS0010	10 rxn (100 µl)
BA01503	100 rxn (1 ml)

Lot no:
Exp. Date:

Storage: -20 °C

Quantity:

The mixture contains 100 units of Atlas HotTaq DNA Polymerase. It is sufficient for 100 reactions of 20 µl or 40 reactions of 50 µl.

Reagents Provided:

- **Atlas HotTaq DNA Polymerase** (0.1 U/µl) in reaction buffer, **4 mM MgCl₂**, **0.4 mM of each dNTPs** (dATP, dCTP, dGTP, dTTP).
- **Water** (nuclease free) - 1 ml.

Description:

Atlas HotTaq 2x PCR Mix is a premixed solution containing everything needed for successful PCR reaction except specific primers and DNA template. The mix includes high-quality recombinant Atlas HotTaq DNA Polymerase, nucleotides and magnesium in a PCR reaction buffer.

For the reaction set-up add the PCR Mix (10 or 25 µl) to the primers, template and water for the total reaction volume of 20 or 50 µl.

To activate the Atlas HotTaq DNA Polymerase it should be incubated at 95 - 97° C for 15 minutes as a first PCR step.

This enzyme allows the PCR setup at ambient temperature without nonspecific annealing and extension.

Quality data:

Activity and stability tested at 20, 30 and 40 cycles of PCR reactions at 95°C. Free of detectable, non-specific nucleases.

Several cycles of freezing/thawing are allowed.

Recommended PCR reaction mix:

Component	Quantity
Atlas HotTaq 2x PCR mix	10 µl (1x)
Primer Forward	0.3 -1 µM
Primer Reverse	0.3 -1 µM
DNA template	1-100 ng/µl
H ₂ O PCR grade	Up to 20 µl
Total	20 µl

Recommended PCR cycles:

Cycle step	Temp.	Time	Cycles
Initial denaturation	95°C	15 min	1
Denaturation	95°C	30-60 s	26-35
Annealing	50-68°C	30-60 s	
Elongation	72°C	1-4 min	
Final elongation	72°C	5-10 min	1

IMPORTANT: Annealing temperature should be 2-6°C lower than the primer melting temperature.

Safety warnings and precautions:

This product is designed for research purposes and *in vitro* use only. According to common laboratory safety practice, it is recommended to wear protective clothing, gloves and safety glasses. Please refer to www.bioatlas.com for Material Safety Data Sheet of the product.

Some applications this product is used in may require a license which is not provided by the purchase of this product. Users should obtain the license if required.