

CappTrio

TRI-VOLUME PIPETTES

As the only pipette manufacturer that offers pipettes with detachable volume controllers, Capp introduces a ready-to-go pipette set which includes three fixed volume knobs.

Highlights:

- Quick volume switch: 3 fixed pipettes in 1!
 - Ideal for routine work
- Customizable Volume Selection
 - Fully Autoclavable



Ordering Information

CappTrio Tri-volume pipettes

Cat. No.	Description	Inaccuracy %	Imprecision %	Colour
C10-T1	Capp Tri-volume 1µl, 5µl and 10µl	2.50 /1.00	1.80 /0.50	○
C100-T2	Capp Tri-volume 20µl, 50µl and 100µl	2.00 /0.70	0.80 /0.20	●
C1000-T3	Capp Tri-volume 200µl, 500µl and 1000µl	1.50 /0.80	0.40 /0.15	●
CXXX-T	Capp Tri-volume Custom Set			