

NEW

CSL-Education Pack

Stock Code: CSL-EDUPACK

New from Cleaver Scientific is a complete low cost DNA electrophoresis Reagent kit.

This kit contains everything required for 20 mini DNA gels in one handy package.

Includes:

- Agarose
- TBE Buffer
- 1 KB Markers
- Runsafe DNA Stain in Sample Loading Dye
- Sample Loading Dye

Contents Include:

- CSL-AG5 Agarose 5 grams x 3
- CSL-MDNA-1KB DNA ladder, 250bp-10KB, 2 x 500ul vial
- CSL-TBEP/2 Powdered Tris-Borate- EDTA Running Buffer - 2 saches (20 litres of 1 x TBE)
- CSL-runSAFE Package 1ml/vial



x 3



x 2

Kit Contains



x 2 Pouches



x 1

Cleaver Scientific Ltd
Unit 4 Triton Park
Brownsover Road
Swift Valley
Rugby
Warwickshire
England CV21 1SG



t +44 (0)1788 565300
f +44 (0)1788 552822
e info@cleaverscientific.com
e sales@cleaverscientific.com
w www.cleaverscientific.com