

Product Data Sheet



Article Nr. LQ-992451

General Information

Description

Bulk, 1000ul Universal Pipette Tips

Raw Material

100% Origin Polypropylene, USP Class VI, free of Heavy metal

Sterilization

Non-sterile

Color

Blue

Compatibility

See Compatibility Chart

Features

Temperature

Approx. -80°C to +121°C

Autoclavability

Autoclave at 121°C/15psi for 20 mins

Shelflife

3 years

Others

Free from Dnase, Rnase, Human DNA, Pyrogen, ATP, Endotoxin

Packaging/Labelling


Packaging

1000/pack 5 packs/ carton

Lot Number

YYMMDD- Year Month Day

Dimension/Drawings



1000µL