

Product Data Sheet

Article Nr. LQ-992450

General Information

Description

Bulk, 200ul Universal Pipette Tips

Raw Material

100% Origin Polypropylene, USP Class VI, free of Heavy metal

Sterilization

Non-sterile

Color

Natural

Compatibility

See Compatibility Chart

Features

Temperature

Approx.-80°C to +121°C

Autoclavability

Autoclave at 121°C/15psi for 20 mins

Shelflife

3 years

Others

Free from Dnase, Rnase, Human DNA, Pyrogen, ATP, Endoxion

Packaging/Labelling


Packaging

1000/pack 20 packs/ carton

Lot Number

YYMMDD- Year Month Day

Dimension/Drawings



200µL