

Serological Pipettes

LabQ Serological Pipettes are manufactured from high-quality virgin polystyrene materials and are suitable for cell culture and microbial culture experiments. The anti-drip tip assures accurate delivery, magnifier stripes and bidirectional graduations make volume reading and size selection easier.



1ml 2ml 5ml 10ml 25ml 50ml 100ml



Printing is crisp, dark, and easy-to-read. Graduation alignment mark allows graduations to be brought into the line of sight quickly. Individually wrapped pipets feature thermoformed paper/plastic wrap to ensure against contamination.

- Six different volumes
- Bidirectional graduations provide choice of ascending and descending scales
- Negative graduations for extra working volume
- Polyolefin plug to help prevent overflow
- Color-coded for quick identification
- Accuracy within $\pm 3\%$ at full volume
- Individually wrapped, paper/plastic
- Nonpyrogenic, Noncytotoxic
- DNase-/RNase-free
- Sterile

Ordering information

Cat.No.	Color Code	Vol. (ml)	Graduations. (ml)	Packaging
LQ-992780	Yellow	1	0.01	50pcs/bag, 10bags/case
LQ-992790	Green	2	0.01	50pcs/bag, 10bags/case
LQ992800	Blue	5	0.1	50pcs/bag, 4bags/case
LQ-992805	Red	10	0.1	50pcs/bag, 4bags/case
LQ-992810	Purple	25	0.2	25pcs/bag, 8bags/case
LQ-992815	Black	50	1	25pcs/bag, 4bags/case